

Rose Bengal test (Brucella test)

- Brucella is a facultative, intracellular, Gram-negative bacteria.
- Brucella species are transmitted from infected animals include goats, sheep, cattle, pigs, and dogs to humans, **The most pathogenic species to humans are:**
- 1- **Brucella melitensis**: Primary host (Goats and Sheep), the most pathogenic in humans.
- 2- **Brucella abortus**: Primary host (Cattle)
- 3- **Brucella suis**: Primary host (Swine-Pigs)
- 4- **Brucella canis**: Primary host (Dogs)

Modes of transmission:

- The consumption of unpasteurized milk and soft cheeses made from the milk of infected animals, especially goats, infected with *B. melitensis* and with occupational exposure of laboratory workers, veterinarians, and butchery workers.

Clinical manifestation:


- ✓ High fever
- ✓ Anorexia
- ✓ Arthralgia (Joint Pain)
- ✓ Tiredness
- ✓ Headache
- ✓ Sweating
- ✓ Lymphadenopathy, Hepatomegaly, and Splenomegaly.

Laboratory Diagnosis (Serological Tests)

Rose Bengal agglutination:

Principle:

Brucella suspended stained with Rose Bengal was used for the direct detection of antibodies in the patient's serum during brucellosis infection.



Kit reagents:


- 1- **Reagent (antigen):** Brucella suspended with Rose Bengal stained.
- 2- **Positive control:** Specific antibodies against Brucella.
- 3- **Negative control:** Non-reactive diluted serum.



A. Slide agglutination test (Qualitative test):

Procedure:

- 1- Bring the reagent and sample at room temperature and mix the reagent vial gently before use.
- 2- Put **1 drop (50 μ l)** of the patient serum on the card.
- 3- Add **1 drop** of Rose Bengal antigen to the serum and mix them.
- 4- Rotate the slide slowly.
- 5- Observe immediately under a light source for any agglutination.
- 6- Report the result.



Read result:

- No agglutination = **Negative (-ve)**
 - Agglutination = **Positive (+ve)**
- ✓ If agglutination appear after **15 seconds = +ve (1:640)**
 - ✓ If agglutination appear after **30 seconds = +ve (1:320)**
 - ✓ If agglutination appear after **1 minute = +ve (1:160)**
 - ✓ If agglutination appear after **1.30 minute = +ve (1:80)**

B. Tube agglutination test (Quantitative test):

- Tube tests are used as quantitative tests to determine the **titer** of antibodies in the samples, by the serial dilution of serum.

For example: If you want to make a serial dilution starting with **(1:20)** and a total volume of **(2 ml)**.

Use this equation: $1/\text{dilution} = \text{amount of solute (serum)} / \text{total volume}$

- ❖ $1/20 = X / 2 \text{ ml}$
- ❖ $X = 2 / 20 = 0.1 \text{ mL (100 } \mu\text{L) amount the serum}$
- ❖ Then add **1.9 ml (1900 $\mu\text{L})$** of diluent (distilled water) to obtain a serial dilution with **(1:20)**.

Procedure:

- 1- Firstly, prepare **(6 or 7)** test tubes.
- 2- To obtain a serial dilution **(1:20)** with total volume **(2 ml)**.
- 3- Add **1.9 ml (1900 $\mu\text{L})$** distilled water + **0.1 mL (100 $\mu\text{L})$** serum to the **first tube only**.
- 4- Add **1 ml (1000 $\mu\text{L})$** distilled water into **other tubes**.
- 5- Take **1 ml** from the **first tube** into the **second tube**, mix, and **repeat up to tube 7**. Then discard the last 1ml from tube 7.
- 6- When serial dilution is finished, add **50 μL** of a reagent (antigen) to each tube.
- 7- A positive reaction is the formation of visible clumps (agglutination) after a period of incubation.

- The **highest dilution** of anti-serum giving agglutination is the **titer** of antibodies.
- The **last tube** that shows an agglutination is called the **end-point** of the test.

