



Washing Cells and Making Red Cell Suspension

- Suspension of washed red cells is needed for all haemagglutination tests.
- The cells must be washed to remove plasma, cell suspensions are used in haemagglutination tests since the ratio of serum to cells affects the sensitivity of most tests – a minimum number of antibodies must bound to RBCs in order to bring about agglutination.

How washing red blood cells:

- 1- Place 0.2 -0.5 ml of blood into the tube (2-3 drops).
- 2- Fill 3/4 of the tube with normal saline (0.9% NaCl).
- 3- Centrifuge for 1-2 minutes until the RBCs are packed.
- 4- Decant the supernatant.



Repeat steps 2-4 at least twice. The last wash should always have a clear supernatant with no signs of hemolysis

Why normal saline (NaCl 0.9%) used in red cell suspension?

- A 0.9% NaCl solution is **isotonic**: when blood cells stay in such a medium, the intracellular and extracellular fluids are in osmotic equilibrium across the cell membrane, and there is no influx or efflux of water.
- In **hypotonic** (e.g. 0.4% NaCl or distilled water), the influx of water occurs that leads to cells swell, membranes disrupted, and hemolysis.



- In **hypertonic** (e.g. 1.8% NaCl), the cells lose their normal shape and shrinks due to the rapid osmotic efflux of water.







Name of the test: Test tube method of blood grouping.

Procedure:

- 1) Collect the blood by vein puncture method.
- 2) Added the blood into the anticoagulant tube (EDTA-tube) and then mixed.
- 3) Added an isotonic solution such as normal saline (N.S) to the mixture.
- 4) Wash the RBCs suspension by centrifugation for 5min. at 3000 Rpm.
- 5) Discard the supernatant and again add normal saline (N.S) to precipitate.
- 6) Mixed the precipitate and repeat steps 4 and 5.
- 7) Re-suspended the precipitate again with saline and standardized it.
- 8) Label three khan tubes with A, B, and Rh (D) and add 2-3 drops of Anti-A, 2-3 drops of Anti-B, and 2-3 drops of Anti-D to tubes A, B, and D. Then, in each of the three tubes, add 2-3 drops of RBC suspension.
- 9) Mix and centrifuge for 1 minute before checking the result.



Result:



↑
Agglutination (positive)

↑
Non-agglutination (negative)