



Rheumatoid Factor Test

- **Rheumatoid factors (RFs)** are **autoantibodies**, usually of the **IgM** class, directed against **human immunoglobulin-G (IgG)**.
- Rheumatoid factors are **elevated** in the blood and synovial fluid of 75% to 85% of people with rheumatoid arthritis (RA).
- It is unclear what triggers the body to produce Rheumatoid factors (RFs).

Principle:

- Testing **latex particles** are coated with **human gamma globulin-G (IgG)** against unknown serums. The presence of a visible agglutination indicates the presence of RF in the samples.



➤ **Reagents:**

- 1- **RF-Latex Reagent:** latex particles coated with human gamma globulin (IgG) in a saline solution.
- 2- **Positive control:** Human serum with an activity about 25 IU/mL.
- 3- **Negative control:** Animal serum with an activity < 5 IU/mL.



Slide agglutination test (Qualitative test):


➤ **Procedure**

- 1- Bring the reagent and sample at room temperature and mix the reagent vial gently before use.
- 2- Put 1 drop (50 μ l) of the patient serum on the card.
- 3- Add 1 drop of RF-Latex Reagent to the serum and mix them.
- 4- Rotate the slide slowly.
- 5- Observe immediately under a light source for any agglutination.
- 6- Report the result.

➤ **Reading:**

Nonreactive: No visible agglutination = **Negative**

Reactive: Any agglutination visible = **Positive**



Negative (-ve)

Positive (+ve)

- **Positive samples** are required to be **re-tested** in **semi-quantitative tests**, which use the same principle and steps as above with only one difference: that the patient serum should be **serially diluted**.

Semi-quantitative Test

Dilute sample with saline as follow:

Tube	Dilution	Composition	RF titer IU/ml
1	1:2	50 µl of serum + 50 µl of saline. Mix	1/16
2	1:4	50 µl from tube 1 + 50 µl of saline. Mix	1/32
3	1:8	50 µl from tube 2 + 50 µl of saline. Mix	1/64
4	1:16	50 µl from tube 3 + 50 µl of saline. Mix	1/128
5	1:32	50 µl from tube 4 + 50 µl of saline. Mix	1/256

Reading:

Titer: is the highest dilution giving agglutination

The result of **RF** (IU/ml) can be obtained by **multiplying** the titer of the dilution by the **minimum** detectable **unit** e.g. titer (1/8).

For example: **Dilution =1:16**

RF concentration= **16 x 8 = 128 IU/mL**